



# Med-Dyne Sleep Study Disposable Electrodes



Leading the Way in Sleep Study Electrode Technology

2775 S. Floyd St.  
Louisville, KY 40209

1-877-MED-DYNE **TOLL-FREE**  
502-429-6759 **FAX**  
sales@med-dyne.com **EMAIL**

[www.med-dyne.com](http://www.med-dyne.com)

# Key Features

- **Available in Cups or Discs** to meet all your EEG needs.
- **Silver/Silver Chloride Electrodes** for preferred economical choice.
- **Color-Coded Lead Wires** made with Industry Standard Safety Sockets, 1.5mm, to ensure compatibility.
- **Customize Color and Length** for your centers' preferences. Choose from Black, Blue, Brown, Green, Orange, Red, White, Yellow, Violet, and Grey in 48", 60", 72", 84", or 96" leadwires.

## Convenient to Use

- 10 Lead sets or individual leads in any length of choice.
- Electrodes work universally with all Colloidal Gels.
- High quality electrodes provide accurate monitoring feedback.
- Save valuable time with single step Morning Disconnect disposal.

## Color-Coded for Easy Connection

- Color-coded lead wires make cable connection easy.
- Strain relief device allows for more connection flexibility.

## Comfortable to Wear

- Latex free.
- Choose Cup or Flat Disc for comfort behind ears.
- Flexible wires feel comfortable against patient's skin.

# Ordering Information



Part #	Description	Quantity
DDJ48	48" Set (BK,WH,GR,RD,BR,BL,OR,YL,VI,GY) EEG Electrode	10 Electrodes /
	Monitoring Disposable Silver/Silver Chloride Discs	1 Set
DDJ60	60" Set (BK,WH,GR,RD,BR,BL,OR,YL,VI,GY) EEG Electrode	10 Electrodes /
	Monitoring Disposable Silver/Silver Chloride Discs	1 Set
DDJ72	72" Set (BK,WH,GR,RD,BR,BL,OR,YL,VI,GY) EEG Electrode	10 Electrodes /
	Monitoring Disposable Silver/Silver Chloride Discs	1 Set
DDJ84	84" Set (BK,WH,GR,RD,BR,BL,OR,YL,VI,GY) EEG Electrode	10 Electrodes /
	Monitoring Disposable Silver/Silver Chloride Discs	1 Set
DDJ96	96" Set (BK,WH,GR,RD,BR,BL,OR,YL,VI,GY) EEG Electrode	10 Electrodes /
	Monitoring Disposable Silver/Silver Chloride Discs	1 Set



Part #	Description	Quantity
DCJ48	48" Set (BK,WH,GR,RD,BR,BL,OR,YL,VI,GY) EEG Electrode	10 Electrodes /
	Monitoring Disposable Silver/Silver Chloride Cups	1 Set
DCJ60	60" Set (BK,WH,GR,RD,BR,BL,OR,YL,VI,GY) EEG Electrode	10 Electrodes /
	Monitoring Disposable Silver/Silver Chloride Cups	1 Set
DCJ72	72" Set (BK,WH,GR,RD,BR,BL,OR,YL,VI,GY) EEG Electrode	10 Electrodes /
	Monitoring Disposable Silver/Silver Chloride Cups	1 Set
DCJ84	84" Set (BK,WH,GR,RD,BR,BL,OR,YL,VI,GY) EEG Electrode	10 Electrodes /
	Monitoring Disposable Silver/Silver Chloride Cups	1 Set
DCJ96	96" Set (BK,WH,GR,RD,BR,BL,OR,YL,VI,GY) EEG Electrode	10 Electrodes /
	Monitoring Disposable Silver/Silver Chloride Cups	1 Set

# Med-Dyne

## Sleep Testing Products




### Gold/Silver Cup Electrodes

- **Available in Gold or Silver Cup Electrodes** to meet your testing preferences.
- **Color-Coded Leadwires** made with Industry Standard Safety Sockets, 1.5mm, to ensure compatibility.
- **Customize Color and Length** for your ECG needs. Choose from Black, Blue, Brown, Green, Orange, Red, White, or Yellow in 48", 60", 72", 96" or 120" leadwires.



### Cloth Monitoring Electrodes


- **Silver/Silver Chloride Snaps.**
- **Breathable Cloth** for long term application.
- **Repositionable Electrodes** without need for added moisture.
- **Hydrogel Electrodes** may be trimmed to fit small babies without loss of adhesive.



Part #	Description	Quantity
DEL3000I-300	7/8" Square Cloth Snap Hydrogel Electrode	300 Electrodes
DEL3000-300	1-1/4" x 7/8" Cloth Snap Hydrogel Electrode	300 Electrodes
DEL3000A-300	1-1/4" Square Cloth Snap Hydrogel Electrode	300 Electrodes

### Foam Monitoring Electrodes


- **Moisture Resistant** to maintain a strong adhesion.
- **Easily Conforms to Skin** with a non-irritating bond.
- **High Quality Foam** paired with a conductive adhesive hydrogel.



Part #	Description	Quantity
DELF3000R-600	1-1/4" Diameter Foam Snap Hydrogel Electrode	600 Electrodes
DELF3000A-600	1-1/4" Square Foam Snap Hydrogel Electrode	600 Electrodes

### Snap Leadwires

- **Compatible with Standard Snap Connectors** for a secure fit to electrodes.
- **Color-Coded Leadwires** made with Industry Standard Safety Sockets, 1.5mm, to ensure compatibility.
- **Customize Color and Length** for your ECG needs. Choose from Black, Blue, Brown, Green, Orange, Red, White, or Yellow in 18", 24", 40", 48", 72", or 120" leadwires.



**Contact your Med-Dyne Representative at 1-877-MED-DYNE to learn more about our product lines.**





# Med-Dyne Value

## Why use Disposable Silver/Silver Chloride Electrodes?

- Significantly reduce your patient's risk of Healthcare-Associated Infections by using new electrodes and lead wires on each patient. 10 Lead sets in any length of choice.
- Americans get an estimated 1.7 million Healthcare-Associated Infections in the U.S. annually, according to the Centers for Disease Control and Prevention (CDC). Each year these infections:
  - Cause about 99,000 deaths.
  - Cost from \$35.7 to \$45 billion in additional healthcare costs annually.

## Why choose Med-Dyne Disposable Silver/Silver Chloride Electrodes?

- Med-Dyne's Disposable Electrodes are designed with Sleep Study centers in mind.
  - Electrodes are the leading economical choice for simplifying Morning Disconnect.
  - Our electrodes work universally with all Colloidal Gels on the market for accurate EEG readings.
  - Med-Dyne has the capability to custom design Electrode sets to fit your needs from color, length, to packaging.
- Med-Dyne always keeps your inventory in stock and ready to ship to reduce your on hand inventory costs.



2775 S. Floyd St.  
Louisville, KY 40209

1-877-MED-DYNE **TOLL-FREE**  
502-429-6759 **FAX**  
sales@med-dyne.com **EMAIL**

[www.med-dyne.com](http://www.med-dyne.com)

Proudly made in the U.S.A.